

NEW NON-INTEGRATIVE MS2-BASED LENTIVIRAL PARTICLES FOR RNA DELIVERY: A SAFE AND EFFICIENT OPPORTUNITY FOR GENE EDITING APPLICATIONS AND THERAPEUTIC PERSPECTIVES

Alexandra Iché¹, Régis Gayon¹, Christine Duthoit¹, Lucille Lamouroux¹, Nicolas Martin¹, Bruno Matéo¹, Jean-Christophe Pagès² and Pascale Bouillé^{1,3}


¹VECTALYS – Canal Biotech II, 3 rue des satellites 31400 Toulouse – France ; ²Inserm U966, Faculté de Médecine, 10 Bd Tonnelier 37000 Tours, France
³FLASHCELL – Canal Biotech II, 3 rue des satellites 31400 Toulouse – France ; Email: tech@vectalis.com

Gene editing using the CRISPR system shows great promises for gene therapy, but it must now face a number of challenges especially for the development of safe and efficient delivery tools for *in vivo*, as well as *ex vivo* gene editing. Indeed, achieving an efficient delivery in hard-to-transfect cells of clinical interest such as T cells, remains challenging and the need for new delivery tools that would allow efficient transfer on most cell types without causing any cell damages is essential for downstream therapeutic applications.


Here, we present an innovative RNA delivery vector, named LentiFlash, based on a new chimeric MS2-based lentiviral platform that allows RNA delivery into the target cells without any genomic scar. This new vector breaks with all existing systems, as the resulting lentiviral particle is able to deliver non-viral coding or non-coding RNA, at high efficiency, into the cytoplasm of any type of target cells.

HIGHLY EFFICIENT DISRUPTION OF PD-1 IN HUMAN PRIMARY T-CELLS


A CRISPR approach using LentiFlash RNA delivery leads to an efficient and dose-dependent knock-out of the PD-1 gene into human activated T lymphocytes, without any genome integration.




LENTIFLASH: A NEW TOOL FOR EFFICIENT CRISPR/CAS9 DELIVERY



LentiFlash allows to achieve more than 80% GFP knock-out in an immortalized monoclonal cell line containing one integrated GFP copy.



THE BENEFITS OF LENTIVIRAL VECTORS WITHOUT ANY GENETIC SCAR



High purification level allows efficient transduction of delicate immune cells without cell phenotype modification.

